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Surecheck®

User information leaflet

As accurate as blood test with 99.9% accuracy on the first day of a missed period. Please read the instructions carefully before use!

For a woman, the prospect of pregnancy is an extremely joyous occasion. Before celebrations, it is advisable to check if you are sure about the pregnancy. Surecheck pregnancy test is developed with special Colloidal Gold Technology, which ensures 99.9% accurate result as early as the first day of the missed period making it as accurate as the blood test. Due to the special Colloidal Gold technology, Surecheck pregnancy test can be used any time of the day and the sleek design of Surecheck makes using the test more hygienic. This makes Surecheck very easy to use. Be assured of accurate result each time, every time.

Store between 4°-30°C

Things to do before you take the Surecheck pregnancy test:

- Check the expiration date on the box. Do not use if expired.
- Have a watch, clock or timer handy.

How to use your Surecheck pregnancy test:

1 Open the protective foil pouch right before use. Remove the test device from the foil package.

2 Urinate on the absorbent tip for 3-5 seconds in order to cover the whole tip.

3 Read the results after 5 minutes.

Your Surecheck box contains:

- A foil pouch with:
 - ✓ Surecheck pregnancy test.
 - ✓ Desiccant (moisture absorbing pouch not required for the test).
- User information leaflet.
- Calendar.

Interpretation of results:

Not pregnant:

Only one color line appears on the control region (C).

Pregnant:

Distinct color lines appear on the control region (C) and test (T) region.

Invalid result:

If no line appears in the control (C) region after 5 minutes, then the result is invalid. It is recommended to repeat the test using a new unopened test kit. The test is rendered invalid if it is removed from the foil pouch for more than 15 minutes before use.

After interpreting your results, do not take any actions without consulting your doctor!

Disregard any changes in the test result after 30 minutes. Surecheck pregnancy test is only for in-vitro diagnostics. Surecheck device should not be reused!

The test can be a wonderful keepsake thanks to the red cap for the absorbent tip and white slider for result window!

FAQs

1. How does Surecheck work?

Surecheck pregnancy test registers the presence of hCG (human Chorionic Gonadotropin), a hormone, in the urine. hCG appears in the urine soon after conception. When urine containing the hCG hormone interacts with the pregnancy test, two lines appear in the result window.

2. When can I use the test?

Surecheck pregnancy test can be used as early as the first day of a missed period. For optimal detection of pregnancy, it is advisable to use the first morning urine sample as it usually contains greater concentration of the hCG hormone. However, due to our special Colloidal Gold Technology, randomly collected urine sample may also be used.

3. What should I do if I have an invalid result (missing control (C) line)? Repeat the process using a new Surecheck pregnancy test.

4. What should I do if the line in the test region (T) is very pale?

Test again using a new Surecheck after 48 hours.

5. What circumstances may affect the result?

- Some drugs may contain hCG hormone. Usage of those products prior to or during the test may affect the result.
- If you use Surecheck 8 weeks following a childbirth or abortion, you may get a false positive result because of the hCG hormone concentration still present in your body. In case you are worried, please contact your doctor to interpret the result.
- Drinking large amounts of liquids may lessen the concentration of the hCG hormone in the urine and may affect the result.

Birth control pills, antibiotics, analgesics (painkillers) and alcohol will not affect the result.

Warning:

Keep out of reach of children. Not for internal use. The test is for self-testing. Test sensitivity: 25 mIU/ml. The results obtained are not to be used for self-treatment. Keep away from direct sunlight, heat and moisture. Surecheck is made for UniComs Switzerland GmbH at Blue Cross Bio-Medical Co. Ltd, Beijiqia industry zone, Beijing, PRC 102209. Store between 4°-30°C.